

CS 550



HIGH DEFINITION SCREEN WITH RED AND GREEN CHARTS



CS 550



CHART SYSTEM

Min : 1.5 m
Max : 6 m



Optometry





HIGH-QUALITY VISUAL ASSESSMENT

- Over 100 kinds of general and special charts for an accurate evaluation of the visual performance
 - Visual acuity: Snellen, Landolt C, symbols for children, numbers and letters including ETDRS
 - Special tests: contrasts, binocular vision, aniseikonia, red/green, astigmatism ...
- Optotypes randomization and reversible tests for mirrored room.
- 2 types of progression available: A type or B type

OPTIMIZED REFRACTION CHAIN

- Can be managed either by remote control or by APH550 automatic phoropter keyboard

OPTIONS

- Floor stand, red/green glass

SCREEN SIZE	19"	TEST CHARTS	Type A (Decimal indication) Type B (Feet and Metric indication)
LUMINOSITY	80 – 320Cd/M ²	Optotypes	Letters, Landolt C, Numbers, Letters, Snellen, Children charts
CONTRASTS	From 100% to 1.5%	Masks and Tests	Astigmatism Binocular vision Coincidence (aniseikonia) Red/Green Heterophoria Stereoscopic vision
PROJECTION DISTANCE	1.5 to 6m		
NUMBER OF CHARTS	100		
PROGRAM SETTINGS	2 specific programmable test sequences		
STANDARD ACCESSORIES	Remote control and mini USB Cable for recharging Wall mounting bracket	DIMENSIONS AND WEIGHT	W 17.2 x H 14.3 x D 2.4 in - 12 lb
		POWER SUPPLY	AV 100-120V/220-240V, 50-60Hz
		 Conformity marking	

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

PRODUCT FUNCTIONS



Screen Size



Projection Distance



Wireless Connection (+IR)



Essilor Instruments USA
8600 W. Catalpa Avenue, Suite 703
Chicago, IL 60656 - USA
Phone: 855.393.4647
Email: info@essilorinstrumentsusa.com
www.essilorinstrumentsusa.com

Your local distributor:

