

SL350 - SL450



ALL THE ESSENTIALS FOR EYE EXAM





SL350 - SL450

SLIT LAMPS

Range of filters



Magnifications



Optometry





WIDE OBSERVATION

- Galilean binocular optical system offering wide field of view and qualitative images.
- Magnification levels easily accessible by drum rotation.



BRIGHT ILLUMINATION

- Tower illumination systems, featuring up to 20° tilting tower illumination.
- Continuous adjustment of slit aperture.
- Available with halogen or LED illumination source.



INTEGRATED BASE

- Easy installation on refraction unit or instrument table thanks to the integrated base and chinrest.

OPTIONS

- Yellow filter, polarised filter, tonometer adaptor.

MAGNIFICATION		SLIT	
Type	Galilean	Slit width	Continuous from 0 to 14mm
Cylinder	12.5x	Slit length	Continuous from 1 to 14mm
Selection	Drum rotation	Aperture	0.2 to 14mm
Field of view	10x (22mm) - 16x (13.5mm) - 25x (8.5mm)	Slit angle	0 to 180°
SL350L/SL350H		Tilt	5, 10, 15, 20°
ILLUMINATION LAMP			
SL450L/SL450H	6x (34.7mm) - 10x (22mm) 16x (13.5mm) - 25x (8.5mm) 40x (5.5mm)	SL350H/SL450H	6V, 20W halogen lamp
		SL350L/SL450L	3v,3W LED lamp
		FILTERS	
BASE		POWER-SUPPLY	
Longitudinal movement	110 mm	220/110V (±10%)	
Lateral movement	110 mm	DIMENSION AND WEIGHT	
Vertical movement	30 mm	450 (W) x 500 (D) x 740 (H) mm 22kg	
Chinrest (vertical)	80 mm	CE Conformity marking	

As improvements are made, these specifications and pictures are not contractually binding and may be changed without prior notice.

Your local distributor:



Essilor Instruments USA
8600 W. Catalpa Avenue, Suite 703
Chicago, IL 60656
Phone: 855.393.4647
Email: info@essilorinstrumentsusa.com
www.EssilorInstrumentsUSA.com