

Cassilo

11580

OPTICAL PRECISION AND HIGH-QUALITY ILLUMINATION ENSURES YOU THE BEST POSSIBLE DIAGNOSIS

## OPH 500 / OPH 700 / RET 500

Range of apertures, graticules and filters



## **OPHTHALMOSCOPES / RETINOSCOPE**

# **OPH 500 - PROFESSIONAL OPHTHALMOSCOPE**

### ULTIMATE CONTROL

- Easier and quicker exams with the positive action maywheel
- Lenses, graticules and filters can all be introduced at the touch of a fingertip

### XENON ILLUMINATION

• Brighter and whiter light. Long lasting bulbs

## COMPREHENSIVE LENS RANGE +29D TO -30D

- Easy examination for all types of ocular structures
- An adaptability for several patient/user prescriptions

#### INCLUDED

• Handle, Lithium battery, case, 1 bulb spare







## OPH 700 - SPECIALIST OPHTHALMOSCOPE THE WORLD'S FINEST DIRECT OPHTHALMOSCOPE • Ergonomic and slim design offers you a confortable fit • The slimline handle is perfectly shaped and balanced HIGH QUALITY ILLUMINATION • The quality optics gives you an undeniable brightness

• Ideal to view through cloudy media and cataracts

## WIDEST RANGE +44D TO -45D

• Wide lens range ensures that you can make optimal diagnosis in any situation

## New Cobalt Blue Filter

• Effective detection and examination of corneal scars and abrasions

#### INCLUDED

• Handle, Lithium battery, case, 1 bulb spare







# **RET 500 - PROFESSIONAL COMBI RETINOSCOPE**

## THE FIRST TRULY BIMODAL RETINOSCOPE

• High performance with a streak retinoscope

### NEUTRALIZATION CHECK

- Express confirmation of neutralization
- Possibility to change the beam from divergent to convergent to check the accuracy of the neutral point

### Two Position Aperture

- Controllable aperture for brightness: 4 mm or 1.7 mm
- Top-quality resolution of the retinal reflex

#### INCLUDED

• Handle, Lithium battery, case, 1 bulb spare







OPHRET 5 KIT KIT OPHTHALMOSCOPE RETINOSCOPE WITH CHARGER

OPH 500 + RET 500 + PSU 500

## OPHRET 7 KIT Premium Kit Ophthalmoscope Retinoscope With Charger

#### OPH 700 + RET 500 + PSU 500



## **PSU 500** LITHIUM-ION DOUBLE CHARGER

- Instrument always ready for use
- Ultimate for long examinations
- LED display for instant indication
- Over two hours of continuous use





#### **OPH 500 Case**



**RET 500 Case** 

CHARACTERISTICS		OPH 500	<b>OPH 700</b>	<b>RET 500</b>
Lens range		+29D to -30D	+44D to -45D	-
Range of apertures			-	4mm to 1.7mm
Range of filters			0 🗖 🗖 🚺 🔵 🚭	
Spare bulbs		OPH 50001	OPH 70001	RET 50001
Brow rest	Long			RET 50002
	Short	<b>_</b>	_	RET 50003
Batteries / Spares		PSU 50004		
Handle sleeves		4 colors: PSU 50005 (blue, orange, green, black)		
Single charger		Mini Lithium charger: PSU 400		
Double charger		Lithium-ion double charger: PSU 500		

## Product apertures, graticules and filters



**Wide angle:** Illuminates the largest area of fundus for the best general diagnosis through a dilated pupil.



**Intermediate:** Easier access through undilated pupils and peripheral examination. Ideal for pediatric examination.



**Macular:** For examination of the macular region of the fundus where a larger beam would create excessive pupillary reaction or patient discomfort.



**Slit:** To determine retinal elevations and depressions; can also be used to assess anterior chamber depth.



**Glaucoma:** Projects a graticule onto the retina to assess the optic cup/disc ratio as an aid to glaucoma diagnosis.

**Fixation Cross:** Projects a graticule onto the retina for assessment of the degree and direction of eccentric fixation.



#### **Red Free and Fixation**

**Cross:** Is for vessel examination with a fixation cross for assessment of eccentric fixation.



**Red Free:** Used to examine blood vessels in fine detail. Blood vessels are silhouetted black against a dark green background.

**Cobalt Blue:** Use with fluorescein dye for examination of corneal scars and abrasions.



Essilor Instruments USA 8600 W. Catalpa Avenue, Suite 703 Chicago, IL 60656 - USA Phone: 855.393.4647 Email: info@essilorinstrumentsusa.com www.essilorinstrumentsusa.com

Your local distributor: