

Small Equipment and Consumables Guide

Specifically designed, developed and tested to enhance lifetime and performance



EDGING PADS



GAM 403

Ø 22 mm • Roll of 1000



GAM 316

Ø 18/14 mm • Roll of 1000



GAM 200

Ø 24 mm • Roll of 1000



GAM 306

Ø 18 mm • Roll of 1000



GAM 207

Lens protectors
Roll of 2000

Compatible Equipment

GAM 222	•			
GAM 218		•		
GAM 216		•		
GAM 202			•	
GAM 203				•

EDGING PLASTIC BLOCKS



GAM 222

Ø 22 mm • Pack of 20



GAM 218

Ø 18/14 mm • Pack of 20



GAM 216

Ø 18/14 mm • Pack of 10



GAM 202

Ø 25 mm • Pack of 20



GAM 203

Ø 18 mm • Pack of 10

Compatible Equipment

MISTER	•	•			
NEKSIA	•	•			
KAPPA/GAMMA			•	•	•
JESS/DELTA	•	•			
SIGMA			•	•	•

CLAMP PADS



M24A33
Ø 22 mm



M1017A
Ø 18/14 mm



N9581D
Ø 25 mm



N9582B
Ø 18 mm



M1012A
Ø 18/14 mm

Compatible Equipment

M24A12	•				
M10A17		•			
M10A09			•		
M10A10				•	
M10A11					•

DRILL/ MILL BITS



FRA210
Drill bit
Ø 1 mm



FRA211
Drill bit
Ø 0.8 mm



FRA308
Drill bit
Ø 0.8 mm



FRA310
Drill bit
Ø 1 mm



M1553A
Drill bit
Ø 0.8 mm



M1554A
Drill bit
Ø 1 mm



FRA212
Mill bit
Ø 1.2 mm



FRA312
Mill bit
Ø 1.2 mm

Compatible Equipment

MR BLUE 2.0			•	•				•
MR BLUE	•	•					•	
MR ORANGE	•	•	•	•				
NEKSIA/KAPPA JESS/DELTA/ SIGMA					•	•		



Wooden Box with 5 Handles and 7 Blades – ECR001

5 nut driver handles ECRM01

	2.2 mm	Square blade	ECRC01
	2.0 mm	One cross blade	ECRC20
	2.0 mm	One star blade	ECRE20
	2.2 mm	One star blade	ECRE22
	2.5 mm	One star blade	ECRE25
	2.2 mm	One hex blade	ECRH22
	1.5 mm	One blade blade	ECRP15



Pocket Wooden Box – TOUR01

			Complete Screwdriver	Spare Blade (x4)
	1.2 mm	Blade screwdriver	TOUP12	TOUL12
	1.6 mm	Blade screwdriver	TOUP16	TOUL16
	1.9 mm	Blade screwdriver	TOUP19	TOUL19
	1.5 mm	Crossblade screwdriver	TOUC15	TOUL15
	2.0 mm	Crossblade screwdriver	TOUC20	TOUL20
	2.5 mm	Crossblade screwdriver	TOUC25	—



Screw Nutdriver Set – TOUA10

	1.6 mm	Blade screwdriver	TOUP36
	1.8 mm	Blade screwdriver	TOUP38
	1.5 mm	Hexa. screwdriver	TOUC35
	1.8 mm	Crossblade screwdriver	TOUC38
	2.02 mm	Hexa. nutdriver	ECRH02
	2.32 mm	Hexa. nutdriver	ECRH32
	2.82 mm	Hexa. nutdriver	ECRH82
	2.32 mm	Star nutdriver	ECRE52
	2.82 mm	Star nutdriver	ECRE82



Mixed Set – TOUA01

	1.2 mm	Blade screwdriver	TOUP10
	1.6 mm	Blade screwdriver	TOUP15
	1.8 mm	Blade screwdriver	TOUP18
	2.3 mm	Blade screwdriver	TOUP23
	1.5 mm	Crossblade screwdriver	TOU006
	2.0 mm	Crossblade screwdriver	TOU007
	2.2 mm	Star nutdriver	ECRE12
	2.5 mm	Star nutdriver	ECRE15
	2.8 mm	Star nutdriver	ECRE18
	2.2 mm	Hexa. nutdriver	ECRH12



MISK08
Tripod



MISK09
Demo lens holder



N82A11
Edger adjusting gauge



L11S44
Pattern holder calibration tool



SPH001
Measuring tools lens Clock



COM003
Large thickness apparatus



L1108P
Touch screen pen



CRA002
Marker bullet tip 1 mm,
opaque ink white lens



99-AOP
Handy All Off Pen
Easily erase progressive
ink markings



99-AOE
Instantly dissolves and removes
ink markings faster and safer than
acetone, in 3 seconds flat – 16 oz.



N82A03
Deblocker plier




JAU002
Lens hole gauge



SET OF 3 REAMERS
1.4 MM – EQU001
1.75 MM – EQU002



EQU003
Set of 2 Needles Reamers
0.5 mm -0.7 mm



EQU010
Set of 6 Reamers
1.4 mm to 2.1 mm



LIMA01
File
Round, thickness: 1.3 mm



LIMA10
Set of 10 Mixed files



TENSISCOPE
TEN070

This tensiscope with tiltable feet makes it possible to identify tensions or stress on mounted lenses.



HAND TENSISCOPE
TEN100

Ultra compact and portable tensiscope for checking stress on mounted lenses

FRAME WARMER

FHE310

Electronic heating regulation for easy frame adjustments



PUPILLOMETER

X81705

Precision measurements with manual control



ULTRASONIC CLEANER

MNU1375/MNU2000

Professional ultrasound cleaning device

MNU2000 is made of stainless steel and holds a 2,000-ml large volume for a busy lab.



OPTEC PLUS SMART VISION SCREENER

C86002/C86004

Enables quick, thorough and repeatable results in pre-testing.

C86004 is the new model with glare and the glare recovery test.



MANUAL LENSMETER

MLM200US

Large ergonomic knobs clearly marked for ease of use.



HIDE-A-BEVEL HAND EDGER

MMWHABHE

Features a diamond wheel with 3 cutting surfaces: 100 grit, 240 grit, and V-Groove.



LENS ENGRAVING DETECTOR

PAL-ID

Easy identification regardless of the lens material



AUTOMATIC LENSMETER

ALM500 / ALM700

Intuitive and practical for rapid and accurate measurement.

ALM700 features measurements of UV transmission and inter-pupillary distance.



TRIAL FRAMES AND TRIAL CASES



Large Trial Lens Set 260 pcs	ESS260
Small Trial Lens Set 167 pcs	ESS167
Trial Frame Metal/Standard	ESS204
Trial Frame Composite/Plastic/Adjustable	ESS304

EXAM LANE ESSENTIALS KIT OPTOESSKIT1

Volk Lens 20D	VLENS20D
Volk Lens 90D	VLENS90D
Horizontal Prism Bar	HORPB
Occluder Pinhole	OCCLUPINM
Flipper Lens +/- 0.25	FLIPLENS025
Original Stereo Fly	SO001
Color Tests (14-plate)	MAG301
Trial Frame Metal/Standard	ESS204
Small Trial Lens Set 167 pcs	ESS167

STEREO AND ISHIHARA TESTS



Original Stereo Fly	SO001
Randot®	SO002
Random Dot "E"	SO003
Butterfly	SO005
Original Randot	SO006
Randot Preschool	SO007
Distance Randot	SO015
Special Randot	SO016
Color Tests (14-plate)	MAG301
Color Tests (24-plate)	MAG302
Color Tests (38-plate)	MAG303

DISPENSING ESSENTIALS KIT DISPESSKIT1

Frame Warmer	FHE310
Pupillometer	X81705
Ultrasonic Cleaner	MNU1375
Lens Engraving Detector	PAL-ID
Manual Lensmeter	MLM200US



Essilor Instruments USA
 8600 W. Catalpa Avenue, Suite 703
 Chicago, IL 60656 - USA
 Phone: 855.393.4647
 Email: info@essilorinstrumentsusa.com
www.essilorinstrumentsusa.com

Your local distributor: