

FOR FAST
AND PRECISE
MEASUREMENT,
JUST LOOK INTO THE
MIRROR





M'EYE FIT MIRROR™

DIGITAL MEASURING DEVICE

M'eye Fit™
MIRROR

A SIMPLE LOOK THROUGH THE MIRROR IS ALL IT TAKES



SIMPLE

Thanks to its compact design and its 30° tilting screen, M'eye Fit Mirror is adapted to all dispensing tables. Plugged to your computer the mirror is immediately operational. The intuitive software will easily guide you through the different steps of measurement.

FAST

No need to raise the camera at eye level, or take photos from multiple angles. With the customer seated at your desk, you achieve precise measurement in two shots. It's simple, fast, and compatible with all types of frames.

PRECISE

M'eye Fit Mirror uses a patented 3D measurement protocol which takes into account the exact angle of the camera and compensates for parallax. A flash secures corneal reflection.

CONNECTION

- Bundled with software
- To be connected to a PC
- Data exchange with shop management software

MEASURES

- Pupillary distance and half-pupillary distance
- Fitting height
- Pantoscopic angle
- Vertex distance

ACCESSORIES

- USB stick
- Storage box
- Frame clip
- Clip stand
- Reading tablet (optional)

DIMENSIONS AND WEIGHT

- 10.6(L) x 7.9(W) x 5.9(H) in
- Weight: 4 lb



ESSILOR
INSTRUMENTS

Essilor Instruments USA
8600 W. Catalpa Avenue, Suite 703 | Chicago, IL 60656 - USA
Phone: 855.393.4647 | Email: info@essilorinstrumentsusa.com
www.essilorinstrumentsusa.com

YOUR CONTACT