



Headrests

Tradition and Innovation – Since 1898, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.



Headrests



No. 10 Stationary

Pillow headrest with sliding concave pad. Adapts to Model 980, 7000 and 6200 chairs as an option. Used primarily as a surgical headrest in the ophthalmic and ENT fields. Standard on the Model 710 chair.



No. 11 Articulating pillow headrest

Designed for the ophthalmic and ENT fields. Standard on the Model 3000 chairs. Used primarily as a diagnostic headrest. Friction lock with adjusting lever.



No. 18 Concave, fully adjustable headrest

Standard headrest on the Model 980, 7000, 6200, and FX-920. Used primarily as a diagnostic headrest in the ophthalmic and ENT fields. Single hand gear lock.



No. 20 Articulating surgical headrest

Adapts to all chairs, except the FX-920, as an option. Used as a surgical headrest in the ophthalmic field. Wrist supports are standard. Friction lock lever.



No. 21 Articulating pillow headrest

Features a sliding concave pad. Adapts to all chairs, except the FX-920, as an option. Used as a diagnostic and surgical headrest for the ophthalmic and ENT markets. Friction lock utilizing lever.



No. 22 Oversized concave, fully adjustable headrest

Adapts to all chairs as an option. Used primarily as a diagnostic headrest for ENT. Single hand gear lock.



No. 27 Special articulating pillow headrest

Features a sliding, removable magnetic pad and a fold-away design for otology. Available on Model 710 only.