

# EYE-RULER™ 2: THE MOST COMPACT MEASUREMENT DEVICE



**QUICK & EASY TO INSTALL**



**WIDE COMPATIBILITY  
WITH IPAD®**



**ULTRA-LIGHT**



**EASE OF USE**




**LATEST EMBEDDED TECHNOLOGIES:  
Bluetooth® low energy**



©Essilor International - September 2019. Essilor® and Eye-Ruler™ 2 are trademarks of Essilor International. ©2019 Apple Inc. All rights reserved. Apple, the Apple logo, iPad® are trademarks of Apple Inc., registered in the U.S. and other countries. The Bluetooth® word mark and logos are registered trademarks owned by the Bluetooth® SIG, Inc. and any use of such marks by Essilor International is under license. Other trademarks and trade names are those of their respective owners.



**ESSILOR  
INSTRUMENTS**



**EYE-RULER™ 2**

**MEASUREMENTS MADE EASY**



Seamless navigation



Personalized lenses



Compact measuring device


**EYE-RULER™ 2**

*For more information,  
please contact your Essilor Instruments representative.*

855.393.4647 | [info@essilorinstrumentsusa.com](mailto:info@essilorinstrumentsusa.com) | [www.essilorinstrumentsusa.com](http://www.essilorinstrumentsusa.com)

Eye-Ruler™ 2 is an all-in-one device enabling personalized measurements with only two photos, as well as frame comparison.

## EYE-RULER™ 2: PERSONALIZED MEASUREMENTS




### FAR VISION MEASUREMENTS

- › Monocular pupillary distance
- › Fitting heights
- › Eye-lens distance
- › Pantoscopic tilt
- › Wrap angle



### NEAR VISION MEASUREMENTS


- › Reading distance



### MEASUREMENT RESULTS

- › User-friendly display with all results at a glance

## EYE-RULER™ 2: FRAME COMPARISON



- › Select and compare up to 4 frames