

Measurements made easy.



A digital measuring application for precise lens centration.




1. Embrace a digital in-store experience

Seamless Navigation

Quickly master fitting parameters and enjoy an easy measurement process thanks to its intuitive interface.

By following the on-screen guidance, the digital display provides reliable results for easy reading and analysis.



For your high-tech practice


The Frame Comparison Module allows the patient to compare and choose their favourite frame by taking pictures of their try-ons to facilitate comparison.



eye-ruler™ 2 is a well-rounded digital measuring application.

With its intuitive navigation, eye-ruler™ 2 is designed to take precise measurements in less than one minute.

It provides practitioners with a complete solution to expertly deliver personalized lenses and ensure a seamless optical experience for the consumer.



2. Benefit from all necessary measurements

For all far vision measurements

Thanks to 3D reconstruction technology, all standard and personalized measurements¹ in less than one minute.

Aim for high precision

Based on efficient and robust hardware and software elements, eye-ruler™ 2 is a compact and handheld solution.

With reliable accessories

- Measuring clip and its correct positioning on the patient's frame
- Ultra-light Bluetooth® add-on device



¹Complete list of all parameters available at the end of this brochure.

3. Raise your game with Near Vision Behavior¹

A unique² lens personalization experience

Exclusive³ to Varilux® XR Series™ lenses, the Near Vision Behavior measurement assesses the patient's visual needs and optimise the shape and position of vision volume.



A 4-in-1 in-store measurement

The natural posture and visual behavior in near vision of the patient is captured by recreating a dynamic reading task of daily life.

Varilux® XR Series™, a tailor-made design from the combination of four parameters for a smooth and complete adaptation to the new lenses.

¹Only available with Varilux® XR Track, XR 4D and XR. Pro.

²Only available to Eye Care Professionals dispensing Essilor® lens innovations compatible with an Essilor Instruments digital device to do the measurements.

³Exclusivity refers to all those products, measurements and technologies covered by patents owned by Essilor®.

Technical specifications

Accessories

eye-ruler2 device



Essilor clip



Personalized measurements

Fitting parameters

Pupillary distances & Half PDs in Far Vision

Fitting heights in Far Vision

ELD (Eye Lens Distance = vertex distance)

Pantoscopic tilt

Wrap angle

Additional parameters

Reading distance

Effective lens diameter

Segment heights

Near Vision PDs

Head cape

Near Vision Behaviour

Hardware specifications of the device

Dimensions	(W) 1.84 x (H) 3.16 x (D) 0.49 in.
Weight	0.13 lbs
Mounting	Magnet
Battery power	2 - 5 days of autonomy
Connection	Bluetooth® low energy
Charger	Magnetic power connection
App compatibility	All iPads with iOS 15 plus



eye-ruler2 Device-Brochure-EN-V2-
Oct2024

ESSILOR INSTRUMENTS USA
8600 W. Catalpa Ave., Ste. 703
Chicago, IL 60656
USA

Made in France.

*As improvements are made, information and data in this product presentation are not contractually binding and may be modified without prior notice. Please check the latest version with your local Essilor Instruments representative.

www.EssilorinstrumentsUSA.com
info@essilorinstrumentsusa.com
855.393.4647

© Essilor International - October 2024

