

A versatile and robust stand and chair
for your clinical needs.



OSC™ 600

Ophthalmic Stand & Chair



1. Proven design

Discover a versatile and robust ophthalmic stand and chair for your clinical needs. The OSC™ 600 integrates all standard functionalities into a compact form, ensuring reliability and ease of use in demanding healthcare environments.



Replaceable upholstery

Offers a practical solution for maintaining cleanliness and comfort in clinical settings.

Robust materials

Constructed with materials designed to ensure durability and reliable performance over time.

Precise movement and positioning

4-axis counterbalanced slit lamp arm provides precise movement and consistent positioning for enhanced operational efficiency.



2. Heavy duty hydraulic lift system

500 lbs patient lift capacity.

Hydraulic lift system ensures controlled and stable adjustments for optimal patient positioning.

Built for longevity engineered with a focus on durability to withstand frequent use in demanding clinical environments.

Manufactured in the USA. Branded with the Essilor Instruments logo.



3. Compact and functional

Non-reclinable minimalist design, optimized for efficient space utilization and practical functionality.

Left or right-handed installation possible, offering flexible configuration to suit any workspace.



Three charging wells providing a convenient option of powering handheld diagnostic devices.

Chair elevation control convenient command on the slit lamp arm.



Technical specifications

Dimensions (chair)

Width	26" / 66cm
Depth	42" / 106cm from footrest
Height	50" / 127cm
Weight	290lbs / 141kg

Environmental conditions in usage (chair)

Temperature	60-104°F / 15-40°C
Humidity	5-95%

Electrical specs (chair)

Power supply voltage	120V
Power supply frequency	60Hz
Amperage	9A
Wall cable length	8ft / 243cm

Chair motorized system (chair)

Type of lift assembly	Hydraulic
Elevation command	Yes - from footswitch and slit lamp arm
Max. patient weight capacity	500lbs / 226kg

Wheelchair compatibility

Compatibility	Yes (optional chair remover required)
---------------	---------------------------------------

Dimensions (stand)

Width	20" / 50cm
Depth	20" / 50cm
Height	86" / 218cm
Weight	291lbs / 132kg

Environmental conditions in usage (stand)

Temperature	60-104°F / 15-40°C
Humidity	5-95%

Electrical specs (stand)

Power supply voltage	120V
Power supply frequency	60Hz
Amperage	10A
Wall cable length	8ft / 243cm

Slit lamp system (stand)

4-axis rotation	Yes
Manual locking mechanism	Yes
Table height range	25.6"-37" / 660-940mm
Max. arm weight capacity	60lbs / 27.2kg

Stand instruments system (stand)

Lamp type	Halogen
Handheld instruments charging wells	3
Phoropter arm height range	10" / 254mm
Electrical sockets for slit lamp on arm	Yes
Electrical sockets for phoropter power supply	Yes
Phoropter cable integration in phoropter arm	Yes

Chair motorized system (stand)

Elevation command control on arm	Yes
----------------------------------	-----

OSC™ 600 is a medical device intended for Optometry. Manufacturer: Haag-Streit USA.

For professional use only, read attentively the instructions for use. As improvements are made, these specifications are not contractually binding and may be modified without prior notice.

OSC 600-Leaflet-US-V1-June2025

OSC™ 600 is a trademark of Essilor International. ©Essilor International - June 2025



NETWORK



(855) 393-4647 | info@essilorinstrumentsusa.com | EssilorinstrumentsUSA.com